

Coverall

Type
5B




Features

Material	60 gsm Microporous film
Fixing	Single / Double stitch / Ultrasonic binding
Cuff	Elastic
Presentation	Standard folding / Medical folding - 40 x 40 cm medical towel optionally wrapped
Sterilization	Non-sterile

- ▶ Optional shoe covers

Standards

- ▶ EN ISO 13982-1
- ▶ EN 14126




Normal



Advanced

Size Chart

/ cm	A	B	C	D
XS	56	156	62	68
S	58	158	64	70
M	60	160	66	72
L	62	162	68	74
XL	64	164	70	76



Packaging

Individual pack

Ordering Code

6041