

# Isolation Gown



## Features

Material	<b>Water-repellent :</b> - 28 gsm hydrophobic PP
	<b>Fluid-resistant :</b> - 35 gsm hydrophobic PP with laminated PE
Fixing	Single / Double stitch
Cuff	Knitted / Elastic
Presentation	Standard folding
Sterilization	Non-sterile

## Standards

- ▶ EN 13795
- ▶ AAMI Level 1, 2, 3

## Colors



Yellow



Blue

## Packaging


- 10 pcs / pack
- 10 packs / case

## Ordering Code


6011

## Size Chart

/cm	A	B	C
XS	106	70	51
S	110	72	58
M	120	74	65
L	123	76	72
XL	128	78	79
XXL	132	80	86



Front



Back

