

Surgical Gown



Features

Material	SMS - 45 gsm hydrophobic PP
Fixing	Single / Double stitch / Ultrasonic binding
Cuff	Knitted / Elastic
Belt	Two ties inside & outside
Presentation	Medical folding - 40 x 40 cm medical towel optionally wrapped
Sterilization	EO

- ▶ Neck-Velcro Closure
- ▶ Four ties system
- maintain sterility of surgeon's backside

Standards

- ▶ EN 13795
- ▶ AAMI Level 1, 2, 3

Advantages

- ▶ Breathable
- ▶ Fluid-repellent
- ▶ Antistatic
- ▶ Low linting

Color



Blue


Packaging

Individual pack


Ordering Code

6021

Layout & Size




Front




Back

Normal



Front



Back

Advanced

Size Chart

/cm	A	B	C
XS	106	70	51
S	110	72	58
M	120	74	65
L	123	76	72
XL	128	78	79
XXL	132	80	86

