



TECHNICAL DATA SHEET

Isolation Gown (6011) AAMI BP:70 Level 1 and 2

Intended Use

Non-surgical Isolation Gown (with belt) intended to protect health care patients and personnel from the transfer of microorganisms, body fluids (Fluid repellent) and particulate material. Single-use and disposable. Latex Free.


Characteristic	Specification			
Description	Gown made of non-woven fabric with back closure, waist and neck ties. Features sleeves closed with an elastic and belt on the right side.			
Material Composition	PE laminated on Polypropylene. Latex-Free			
Size	Regular (120cm x 140cm)			
Color	Blue			
Finishing	Level 1		Level 2	
	Stitch sewing		Ultrasonic sealing	

Description	Dimension/Size(cm)
Length(A)	120
Chest Circumference(B)	140
Sleeve Length (C)	59
Cuff(D)	5
Muffs(E)	28
Shoulder(F)	25
Collar(G)	22
1、 Lacing of the neckline 40±3cm 2、 Waistbelt 170±10cm centered 3、 Cuff(terylene)8x5cm	




TECHNICAL DATA SHEET

Isolation Gown (6011)



Front



Back

Ordering Information

Description	EAN-13	Pcs/Bag	Boxes/Case	Each/Case	Carton Size	Netweight
Isolation Gown, Regular Size	4897030531095	10	10	100	59*26*43cm	N.W.:9.1KG G.W.:10.1KG



TECHNICAL DATA SHEET

Isolation Gown (6011)

Visual



Storage Condition: Store under ambient temperature and dry condition
Country of Origin: China



Hong Kong Office:
SKYPRO MEDICAL SUPPLIES COMPANY
Room 301-303, Block C, Phase 2,
Tsing Yi Industrial Centre, 1-33 Cheung Tat Road
Tsing Yi, New Territories, Hong Kong

Tel: +852 - 2711 0882
E-mail: info@1masks.com



Factory :
SPRO MEDICAL PRODUCTS(XIAMEN)CO.,LTD
No.139, TongAn Garden, TongAn Industry Area,
Xiamen City, Fujian Province, 361100, P.R.China

Tel: 86-0592-3697888
E-mail: info@1masks.com