

# Technical Specification Sheet

## Surgical Mask Type IIR

Type IIR medical face masks should only be used for patients and other persons to reduce the risk of spread of infections particularly in epidemic or pandemic situations.

Surgical Face Mask Type IIR		
Characteristic	European Standard	Acceptable Limits
Bacterial Filtration Efficiency(BFE)	EN14683 IIR	BFE @ 3.0 $\mu$ m large Bacteria $\geq$ 98%
Differential Pressure (Delta P)		$\Delta$ P < 6.0 mm H <sub>2</sub> O/cm <sup>2</sup>
Synthetic Blood Penetration		120mmHg
Microbial cleanliness		< 30 cfu/g
Flammability		Class I flame spread

Description	Specification		
Flat pleats rectangular facemasks with a shapeable nose-piece and two earloops present, one on each side, in order to hold mask in place.	Mask Body	Outer Material	Spunbond (SBPP)
		Filter Layer	Melt Blown (MBPP)
		Inner Material	Spunbond (SBPP)
	Nose-piece	PE or Metal wire covered by PE	
	Earloop	Polyester and Spandex	
	Body Size	Length	175mm $\pm$ 3mm
		Width	95mm $\pm$ 2mm
	Length of Earloop	Regular	145mm $\pm$ 3mm
Length Nose-piece	125mm $\pm$ 2mm		

Packaging Information	Color	Origin
50 Masks/Box, 40 Boxes/Case, 2000 Masks/Case	Blue	China

# Technical Specification Sheet

## Surgical Mask Type IIR

### Visual



### Ordering Information

Description	EAN-13	Pcs/Box	Boxes/Case	Each/Case	Carton Size	Neiweight
Surgical Masks Type II R	4897030530845	50	20	2000	52.5*38* 42.5cm	N.W:7KG G.W:9.8KG



Hong Kong Office:  
 SKYPRO MEDICAL SUPPLIES COMPANY  
 Room301-303, Block C, Phase2,  
 Tsing Yi Industrial Centre, 1-33 Cheung Tat Road  
 Tsing Yi, New Territories, Hong Kong

Tel: +852 - 2711 0882  
 E-mail: info@1masks.com



Factory :  
 SPRO MEDICAL PRODUCTS(XIAMEN)CO.,LTD  
 No.139, TongAn Garden, TongAn Industry Area,  
 Xiamen City, Fujian Province, 361100, P.R.China

Tel: 86-0592-3697888  
 E-mail: info@1masks.com