



SkyPro

# **PRODUCT CATALOGUE**

MEDICAL PRODUCTS



Quality First

Customer First

# ABOUT SKYPRO

SKYPRO Medical Supplies Company is a Hong Kong-based enterprise founded in 1998. We have been a supplier of medical consumables for Hong Kong Hospital Authority since 2003 and established a manufacturing base in Xiamen, Fujian Province of China in 2006. Our business covers the R&D, design, manufacture, marketing, and sales of high-quality single-use medical masks, respirators, protective clothing, and drapes. Our products have been sold to more than 40 countries and regions, occupying a large market share.



## Quality First! Customer First!

**"We are committed to provide consistent, high quality medical products for the community and engage with our customers for improvement."**


**"Quality First!"** – Our company has implemented strict, effective, systematic manufacturing and management systems, and obtained ISO13485, ISO9001 and CE certificates, national industrial production licenses, and GMP qualifications for standardized production. In addition, our R&D teams in Hong Kong and Xiamen pay tremendous efforts in the material selection and engineering to ensure our products always at the top standard.

**"Customer First!"** – Our company provides a wide range of customer services to meet our customers' needs, and set up a three-dimensional sales and marketing system, "Original Equipment Manufacturer (OEM) / Original Design Manufacturer (ODM) + E-Commerce + Flagship Store".

## MANUFACTURING PLANT



Gown manufacturing section




Dressing manufacturing section




Mask manufacturing section



Mask quality control section



Product output section



Sterilization section

## CERTIFICATES





# CONTENT

MEDICAL  
MASK

\_\_\_\_\_ 4

RESPIRATOR

\_\_\_\_\_ 8

GOWN

\_\_\_\_\_ 15

DAILY MASK

\_\_\_\_\_ 9

COVERAGE



\_\_\_\_\_ 20



# Procedure Mask




## Features

EN 14683 Type*	I : BFE $\geq$ 95% II : BFE $\geq$ 98%
Material	PP
Inner Layer	Hydrophobic  Hydrophilic 
Ply	3
Fold	3 / 4
Tie Type	Earloops





\* Tested by Nelson Labs

## Colors



▶ Customized color or print is also available

## Ordering

	Code	Fold	Inner Layer
Type I	1131	3	
	1141	4	
	1132	3	
	1142	4	
Type II	1231	3	
	1241	4	
	1232	3	
	1242	4	

## Sizes

Adult	175 x 95 mm
Child	145 x 95 mm
Kid	120 x 85 mm

## Packaging



5	pcs / zip bag
50	pcs / box
40	boxes / case

▶ Optional individual pack

# Activated Carbon Mask



## Features

EN 14683 Type*	I : BFE $\geq$ 95% II : BFE $\geq$ 98%
Material	PP
Inner Layer	Hydrophobic  Hydrophilic 
Ply	4
Fold	3 / 4
Tie Type	Earloops





\* Tested by Nelson Labs

## Color



Carbon Black

## Ordering

	Code	Fold	Inner Layer
Type I	1133	3	
	1143	4	
	1134	3	
	1144	4	
Type II	1233	3	
	1243	4	
	1234	3	
	1244	4	

## Size

175 x 95 mm

## Packaging

50 pcs / box

40 boxes / case

► Optional individual pack

# Surgical Mask



## Features

Material	PP
Inner Layer	Hydrophobic
Ply	3
Fold	3
Tie Type	Ties / Earloops

## Standards

- ▶ EN 14683 Type I, II, IR, IIR
- ▶ ASTM Level 1, 2, 3

## Colors



Blue



Green



White

## Size

175 x 95 mm

## Packaging

50 pcs / box  
40 boxes / case

## Ordering Code

SP02

# Dental Mask - with Eye Shield

## ASTM Level


1

2

3



## Features

	BFE	PFE	Splash Resistance
1 :	≥ 95 %	≥ 95 %	> 80 mmHg
ASTM Level 2 :	≥ 98 %	≥ 98 %	>120 mmHg
3 :	≥ 98 %	≥ 98 %	>160 mmHg
Material	PP		
Inner Layer	Hydrophobic 		
Ply	3		
Fold	3		
Tie Type	Ties / Earloops		

## Colors



Blue



Green



White

## Size


175 x 95 mm

## Packaging

25 pcs / box

40 boxes / case

## Ordering

	Code	Inner Layer
1	2131E	
ASTM Level 2	2231E	
3	2331E	



# Respirator

KN95

FFP2



## Standards

- ▶ GB 2626
- ▶ EN 149

## Features

Material	PP
Inner Layer	Hydrophobic / Hydrophilic
Ply	3 / 4
Tie Type	Earloops / Straps
Valve	With / Without
Nose Bridge	Enclosed / Exposed
Edge Curve	Sharp angle / Flat angle
Sterilization	Sterile / Non-sterile

\* Tested by Nelson Labs

## Colors



White

Navy Blue

▶ Customized print is also available

## Size

160 x 110 mm

## Packaging

GB 2626	40	pcs / box
	30	boxes / case
EN 149	50	pcs / box
	10	boxes / case

## Ordering Code

D918	
GL001	(With Valve)
GL001A	(Without Valve)

# Catering Mask



## Features

Material	PET
Tie Type	Earloops

## Advantage

- ▶ Prevent food contamination

## Color

- ▶ Transparent with white frame

## Size

200 x 140 mm

## Packaging


Individual pack

## Ordering Code

3001





# Foldable Mask - Adjustable Earloops



## Features

EN 14683 Type*	I : BFE $\geq$ 95% II : BFE $\geq$ 98%
----------------	---

Material	PP
----------	----

Inner Layer	Hydrophobic  Hydrophilic 
-------------	--

Ply	3
-----	---

Tie Type	Earloops
----------	----------

\* Tested by Nelson Labs

## Advantage

- ▶ Two wearing methods
  - Earloops
  - Headbands

## Color



White





## Size

185 x 80 mm

## Packaging

20	pcs / box
40	boxes / case



## Ordering

	Code	Inner Layer
Type I	1121	
	1122	
Type II	1221	
	1222	

# Foldable Mask - Stretchable Non-woven Earloops



## Features

EN 14683 Type*	I : BFE $\geq$ 95% II : BFE $\geq$ 98%
Material	PP
Inner Layer	Hydrophobic  Hydrophilic 
Ply	3
Tie Type	Earloops

\* Tested by Nelson Labs





## Advantage

- ▶ Comfortable earloops

## Color

- ▶ Customized color or print

## Ordering

	Code	Inner Layer
Type I	1111	
	1112	
Type II	1211	
	1212	

## Size

135 x 90 mm

## Packaging

20	pcs / box
40	boxes / case

# Reusable Mask - Cotton



## Features

Material	Cotton
Ply	2
Tie Type	Earloops

- ▶ With replaceable filter

## Advantage

- ▶ Washable

## Color

- ▶ Customized color or print



## Size

200 x 140 mm

## Packaging


Individual pack

## Ordering Code

3012



# Reusable Mask - PU



## Features

Material	PU
Ply	1
Tie Type	Earloops



## Advantage

- ▶ Washable

## Color

- ▶ Customized color



## Sizes

Adult	160 x 120 mm
Child	145 x 114 mm

## Packaging

- 3 pcs / zip bag
- With individual pack



## Ordering Code

3014

# Sport Mask



## Size

280 x 160 mm

## Packaging

Individual pack

## Ordering Code

3011

## Features

Material	Neoprene fabric
Ply	1
Tie Type	Velcro headband

► With replacable filter

## Color

► Customized print

# Isolation Gown



## Features

Material	<b>Water-repellent :</b> - 28 gsm hydrophobic PP
	<b>Fluid-resistant :</b> - 35 gsm hydrophobic PP with laminated PE
Fixing	Single / Double stitch
Cuff	Knitted / Elastic
Presentation	Standard folding
Sterilization	Non-sterile

## Standards

- ▶ EN 13795
- ▶ AAMI Level 1, 2, 3

## Colors



Yellow



Blue

## Packaging

- 5 pcs / pack
- 10 packs / case

## Ordering Code

6011

# Surgical Gown



## Features

Material	SMS - 45 gsm hydrophobic PP
Fixing	Single / Double stitch / Ultrasonic binding
Cuff	Knitted / Elastic
Belt	Two ties inside & outside
Presentation	Medical folding - 40 x 40 cm medical towel optionally wrapped
Sterilization	EO

- ▶ Neck-Velcro Closure
- ▶ Four ties system  
- maintain sterility of surgeon's backside

## Standards

- ▶ EN 13795
- ▶ AAMI Level 1, 2, 3

## Advantages

- ▶ Breathable
- ▶ Fluid-repellent
- ▶ Antistatic
- ▶ Low linting

## Color



Blue

## Packaging


Individual pack

## Ordering Code


6021

# Gown Layout & Size

## Isolation Gown




Front




Back

## Surgical Gown




Front




Back

Normal



Front




Back

Advanced

## Size Chart

/cm	A	B	C
XS	106	70	51
S	110	72	58
M	120	74	65
L	123	76	72
XL	128	78	79
XXL	132	80	86





# Scrub Suit



## Features

Material	72 gsm Spunlace
Fixing	Single / Double stitch
Cuff	Knitted / Elastic
Presentation	Standard folding
Sterilization	EO

## Color



Blue

## Packaging


Individual pack

## Ordering Code

6031

## Size Chart

/cm	A	B	C	D	E	F	G
XS	66	49	19	99	53	37	31
S	69	52	20	104	57	41	33
M	72	54	20	108	60	42	35
L	74	55	21	111	61	45	36
XL	77	55	21	116	62	45	37
XXL	80	58	22	118	63	48	38



# Coverall





## Features

Material	60 gsm Microporous film
Fixing	Single / Double stitch / Ultrasonic binding
Cuff	Elastic
Presentation	Standard folding / Medical folding - 40 x 40 cm medical towel optionally wrapped
Sterilization	Non-sterile

► Optional shoe covers

## Standard


► Type 5 - EN ISO 13982-1



GOWN

## Size Chart

/ cm	A	B	C	D
XS	56	156	62	68
S	58	158	64	70
M	60	160	66	72
L	62	162	68	74
XL	64	164	70	76



## Packaging

Individual pack

## Ordering Code

6041

# Coverage

## ► Surgical Cap



### Features

Material	25 gsm PP
Fixing	Single stitch
Closure	Ties / Elastic Band / Velcro
Size	One-size-fits-all

### Ordering Code

7011

## ► Strip Bouffant Cap



Material	10 gsm PP
Fixing	Double stitch
Size	21 inches

7021

## ► Round Cap



Material	20 gsm PP
Fixing	Single stitch
Size	21 / 23 / 25 inches or customized

7031

### Packaging

100 pcs / pack  
10 packs / case

### Colors



White



Blue

# Coverage

## ▶ Sleeve Covers



Advanced

PE  
laminated

Material

30 gsm PP

Colors



Blue



White

Ordering  
Code


7041

## ▶ Boot Covers



Advanced

Non - skid



30 gsm  
Microporous film



White

7051

## ▶ Shoe Covers



Advanced

Non - skid



PE  
laminated

30 gsm PP



Blue



White

7061

Packaging

100 pcs / pack

10 packs / case

# Coverage

► Bed Sheet



Advanced

PE  
laminated

Material

30 gsm PP

Colors



Blue



White

Ordering  
Code

7071

► Pillow Case



Advanced

PE  
laminated

30 gsm PP



Blue



White

7081

Packaging

Customized pack





# Medical Mask

## Billing Information

Name : \_\_\_\_\_ Address : \_\_\_\_\_  
Company : \_\_\_\_\_  
Phone No. : \_\_\_\_\_ City : \_\_\_\_\_  
Email : \_\_\_\_\_ State : \_\_\_\_\_ Zip Code : \_\_\_\_\_

Item Code : \_\_\_\_\_

Quantity : \_\_\_\_\_

## Size

- ☐ 175 x 95 ± 2 mm  
☐ 145 x 95 ± 2 mm  
☐ 120 x 85 ± 2 mm

## Color

- ☐ White ☐ Black  
☐ Green ☐ Blue  
☐ Yellow ☐ Pink

## Graphic Printing

- ☐ Yes  
Remarks : \_\_\_\_\_  
☐ No

## Embossing Logo

- ☐ Yes - Please mark the logo position  
☐ No



## Nose Bridge

- ☐ PE-coated single wire  
☐ PE-coated double wires  
☐ PE  
☐ Others : \_\_\_\_\_

## Tie Type

- ☐ Ties  
☐ Earloops { ☐ Flat  
                                  ☐ Round

Remarks : \_\_\_\_\_

## Packaging

- ☐ Individual polybag  
☐ \_\_\_\_\_ pcs / zip bag  
☐ \_\_\_\_\_ pcs / box

## Fabric Composition

- ☐ Outer-layer PP : \_\_\_\_\_ gsm  
☐ Middle-layer MB : \_\_\_\_\_ gsm  
☐ Inner-layer\* PP / ES : \_\_\_\_\_ gsm  
\*Hydrophilic (A) / Hydrophobic (R)

## Special Requirements

Signature : \_\_\_\_\_

Date : \_\_\_\_\_



# Respirator

## Billing Information

Name : \_\_\_\_\_ Address : \_\_\_\_\_  
Company : \_\_\_\_\_  
Phone No. : \_\_\_\_\_ City : \_\_\_\_\_  
Email : \_\_\_\_\_ State : \_\_\_\_\_ Zip Code : \_\_\_\_\_

Item Code : \_\_\_\_\_

Quantity : \_\_\_\_\_

## Color

- ☐ White  
☐ Black  
☐ Navy Blue  
☐ Others : \_\_\_\_\_

## Embossing Logo

- ☐ Yes — Please mark the logo position  
☐ No



## Nose Bridge

- ☐ PE-coated double wires  
☐ PE  
☐ Others : \_\_\_\_\_

## Graphic Printing

- ☐ Yes { ☐ Whole area  
                  ☐ Selected position  
☐ No

## Packaging

- ☐ Individual polybag  
☐ \_\_\_\_\_ pcs / zip bag  
☐ \_\_\_\_\_ pcs / box

## Tie Type

- ☐ Headbands  
☐ Flat earloops

Remarks : \_\_\_\_\_

## Fabric Composition

- ☐ Outer-layer PP : \_\_\_\_\_ gsm  
☐ Middle-layer MB : \_\_\_\_\_ gsm  
☐ Inner-layer\* PP / ES : \_\_\_\_\_ gsm

\*Hydrophilic (A) / Hydrophobic (R)

Remarks : \_\_\_\_\_

## Special Requirements

Signature : \_\_\_\_\_

Date : \_\_\_\_\_



info@1masks.com



Hong Kong



+852 2711 0882



Room 301-303, 3/F, Block C,  
Phase 2 Tsing Yi Industrial Center,  
1-33 Cheung Tat Road, Tsing Yi



China



+86 592 319 7888



No.139 Factory Building,  
TongAn Garden, TongAn Industrial Area,  
TongAn, Xiamen

[www.1masks.com](http://www.1masks.com)

SKYPRO Medical Products Co., Ltd.